

Small HCPs in a 12 kDa Protein Drug Analyzed by GeLC-MS/MS

Rikke Raaen Lund¹, Marie Grimstrup¹, Maheen Sayeed², Karin Abarca Heidemann², Ejvind Mørtz¹
¹Alphalyse A/S, Odense, Denmark, ²Rockland Immunochemicals Inc., Limerick, PA, USA



Introduction

The content of low molecular weight host cell proteins (HCPs) in purified protein drugs is often difficult to evaluate, due to their low immunogenicity and poor ability to be visualized in gel-based total protein stains. The proteome of commonly used expression organisms, such as *E. coli* and Chinese Hamster cells, contains 30-40% proteins with a molecular weight below 20kDa, and these are easily missed in both in gel separations, Western blots and ELISA quantitation of the total HCP-content. To provide unbiased analysis of small as well as larger HCPs, we introduce the use of a mass spectrometry-based orthogonal method, well known from proteomics, called GeLC-MS/MS.

Here, we analyze an in process protein drug, a 12 kDa protein produced in *E. coli*, as well as the corresponding null cell lysate, using 1D-PAGE and nano-flow LC-MS/MS (GeLC-MS/MS) to achieve high coverage of small HCPs.

1D-PAGE

Proteins were separated by 1D-PAGE and stained for isolation and visualization of the high concentration protein drug in separate gel fractions. Most proteins, 78% in the null cell lysate and 82% in the in process sample were only identified in a single gel fraction (light green and blue bars, fig. 1 and 2, respectively). Further, the fractionation leads to a high number of protein identifications in the fractions without the protein drug (fig 2).

Protein separation by 1D-PAGE prior to LC-MS/MS leads to identification of a high percentage of small HCPs: 46% of the HCPs identified in the null cell lysate and 37% in the in process sample were smaller than 20kDa (fig. 3, table 1 and 2). Further, the HCPs that were identified covered the entire pI range and 99% of the molecular weight range of the total *E. coli* proteome (fig. 3). Examples of small HCPs is given in table 3 and 4.

The sensitivity of the method was estimated by parallel analysis of the purified protein drug with two spiked-in standards, these standards were both identified at 50ppm.

Nano Flow LC-MS/MS

The proteins contained in the gel fractions were digested enzymatically into peptides using trypsin. Separation of the complex peptide mixture was achieved by nanolitre flow HPLC using a CSH (charged surface hybrid, Waters) column. This material has a high loading capacity and excellent peak shape in formic acid mobile phases. This enables highly sensitivity and accurate MS detection of low level HCPs using a qTOF mass spectrometer (Bruker Corp. fig. 4 and 5).

GeLC-MS/MS advantages

- 1D-PAGE requires very little sample preparation and has a wide mass and pI range, leading to unbiased sample analysis
- Protein digest combined with nano flow LC-MS/MS provides sensitive and accurate measurement of multiple peptides from each protein as well as sensitive and accurate measurement of peptide fragment masses

HCP Identification Criteria

HCPs are identified by comparing the mass data to the theoretical masses of the host cell proteome. Peptides are assigned a score according to how well they match and these are summarized to a protein score. Since small proteins have fewer peptides they also have lower proteins scores even at the same molar level as larger proteins.

We have used conservative inclusion criteria that accounts for this: 1% false discovery rate, a minimum score of 25 for all peptides, and a two peptide minimum for proteins larger than 20kDa. These criteria allows highly confident HCP identification with a low false positive rate.

Table 3: Top 30 small proteins in the null cell lysate

HCP no.	Accession no.	Protein Name	Mass	pI	Score
22	tr C6CEG6	Ribosomal protein L9	15759	6.17	3350
29	tr C6CEG3	Ribosomal protein S7	17593	10.3	2701
31	tr C6CEG9	6,7-dimethyl-8-ribitylmuramoyl synthase	16147	5.15	2656
42	tr C6CEG0	PTS system, glucose subfamily, IIA subunit	18240	4.73	2336

Table 4: Top 30 small HCPs in the in process sample

HCP no.	Accession no.	Protein Name	Mass	pI	Score
60	tr C6CEE3	Ferric uptake regulator, Fur family	17012	5.68	3455
61	tr C6CEG3	Protoporphyrinogen IX synthase	19430	5.23	3385
62	tr C6CEG5	Methionine-S-sulfide reductase	15789	5.58	3381
70	tr C6CEG1	30S ribosomal protein L10	17757	9.04	3367
72	tr C6CEG2	Ribosomal protein S7	17593	10.3	3082
73	tr C6CEG7	Glutaredoxin	13042	4.75	3066
75	tr C6CEG9	Ribosomal protein L29	7269	9.38	3025
78	tr C6CEG9	30S ribosomal protein S5	17534	10.23	2991
79	tr C6CEG8	30S ribosomal protein L11	14923	9.64	2999
82	tr C6CEA7	DNA protection during starvation protein	18684	5.72	2990
83	tr C6CEG0	Regulator of sigma D	18288	5.65	2971
86	tr C6CEG0	Molybdenum cofactor sulfur carrier subunit	8724	4.38	2970
88	tr C6CEA6	Ferritin	19468	4.77	2932
96	tr C6CEA8	Iron-sulfur cluster assembly scaffold protein IscU	14011	4.82	3011
99	tr C6CEG7	30S ribosomal protein L31	8094	9.46	2711
105	tr C6CEG3	30S ribosomal protein S6	15163	5.25	2621
107	tr C6CEG8	ATP synthase F1, delta subunit	19434	4.94	2541
108	tr C6CEG7	Uncharacterized protein GN+ECRD_4172	10323	5.18	2321
110	tr C6CEG8	Flavodoxin	15896	4.21	2241
112	tr C6CEG4	Ribosomal protein L18	12762	10.41	2221
95	tr C6CEG0	DNA-binding protein	15587	5.43	2021
96	tr C6CEG1	30S ribosomal protein L6	18949	9.71	2011
97	tr C6CEK7	UspA domain protein	15925	6.03	2001
102	tr C6CEH0	Acyl carrier protein	8634	3.98	1901
103	tr C6CEA8	Iron-sulfur cluster assembly accessory protein	12264	4.11	1891
105	tr C6CEG5	30S ribosomal protein L22	12210	10.23	1851
111	tr C6CEG9	Ribosomal protein S11	13950	11.19	1581
112	tr C6CEG0	PTS system, glucose subfamily, IIA subunit	18240	4.73	1531
113	tr C6CEG5	Acetyl-CoA carboxylase, biotin carboxyl carrier protein	16733	4.66	1521
115	tr C6CEK4	Thioredoxin	15887	5	1471

Figure 4: MS/MS spectrum of the a peptide from Ferric Uptake Regulator (in process sample)




Figure 5: CSH column and qTOF Mass spectrometer



Figure 3: Molecular weight and pI of the HCPs identified by GeLC-MS/MS and the E.coli proteome




Table 1: Distribution of high and low molecular weight HCPs in the null cell lysate

Protein size	HCPs identified
HMW (≥20kDa)	297
LMW (<20kDa)	256
Total	553

Table 2: Distribution of high and low molecular weight HCPs in the in process sample

Protein size	HCPs identified
HMW (≥20kDa)	96
LMW (<20kDa)	56
Total	152

Figure 1: GeLC-MS/MS of null cell lysate





Figure 2: GeLC-MS/MS of the in process sample



HCPs Characteristics and Individual Risk Assessment

Identification of the protein names and database accession numbers enables an in-depth analysis of each individual HCP present in the drug sample. Important features for a drug risk assessment include:

- Immunological properties and presence of human B- and T-cell epitopes
- Homology to human proteins with important biological function
- Homology to the drug protein
- Enzymatic activity to modify or cleave the drug product constituents
- Hormone or hormone-like activity

Known information can be found at <http://www.uniprot.org/> including:

- Molecular function, biological process, ligand binding and cellular component
- Post translational modifications and processing
- Expression, interaction and structure
- Protein family and domains
- Sequence variations
- Publications

Prediction of protein features can be investigated at <http://www.cbs.dtu.dk/services/>, including (fig. 6):

- Immunological features
- Post translational modifications
- Protein structure and function

Figure 6: Examples of services available for prediction of protein features

Immunological features	Post-translational modifications of proteins
Agnitype	Dalysite
Residue-level epitope mapping of antigens based on peptide microarray data	C-glycosylation sites (trained on Dicytostelium discoideum proteins)
Near B-cell epitopes	Discontinuous B-cell epitopes
HLA binding	Hydrophobicity
Antibody recognition	Hydrophobicity
N-glycanase	Identification of Rossmann folds and prediction of FAD, NAD and NADP specificity
Proteinase K resistance	Identification of Rossmann folds and prediction of FAD, NAD and NADP specificity
Proteinase K resistance	Identification of Rossmann folds and prediction of FAD, NAD and NADP specificity
Proteinase K resistance	Identification of Rossmann folds and prediction of FAD, NAD and NADP specificity
Proteinase K resistance	Identification of Rossmann folds and prediction of FAD, NAD and NADP specificity
Proteinase K resistance	Identification of Rossmann folds and prediction of FAD, NAD and NADP specificity

Conclusion

The HCPs that were identified covered 99% of the entire *E. coli* proteome in terms of molecular weight and pI. This shows that the developed GeLC-MS/MS method has no inherent limitations with respect to pI or molecular weight for HCP identification. The obtained protein identity enables an in-depth analysis of each individual HCP and a more detailed risk assessment.